


CRBasic
SI-111 1969 Custom Coefficients

Calibration Date: 2-Jun-2008

	C2	C1	C0
m	56229.3	9.16E+06	1.49E+09
b	29536.7	-8.52E+05	-2.00E+07

Use these coefficients under "Declare Constants" section in the example Campbell Scientific datalogger program



Edlog
SI-111 1969 Custom Coefficients

Calibration Date: 2-Jun-2008

	C2	C1	C0
m(SB)	0.56230	91.62252	14897.54898
b(SB)	0.29537	-8.51635	-200.07200

Use these coefficients in instructions 11&13 in the example Campbell Scientific datalogger program