

QSR-240 Quantum Scalar Reference Sensor

The QSR-240 Hemispherical Quantum Scalar Reference Sensor is designed for monitoring total incident Photosynthetically Available Radiation (PAR; 400-700 nm) from the sun and sky. The rugged QSR-240 features a patented solid Teflon® spherical collector with hemispherical field-of-view cutoff plate, ensuring uniform directional response over 2 π steradians. An aluminumencased optical light pipe funnels flux from the collector to a silicon photodetector that has a flat quantum response over PAR. The QSR-240 sensor shares the scalar directional response characteristics of the Biospherical Instruments QSP-200 series scalar sensors

QSR-240 applications include monitoring solar radiation in energy studies and aquatic productivity studies, and as a surface irradiance reference sensor for the QSP-200 series profiling sensors or PNF-300 Profiling Natural Fluorometer.

Options for the QSR-240 include QSC-240 weather-resistant cable, DAS-186 automated Data Acquisition System and power supply, or QSP-170BD digital display and power supply.



Key Features

- Designed to monitor total incident PAR (400-700 nm) from the sun and sky
- 1.9 cm (3/4") diameter Teflon® spherical irradiance collector with hemispherical field-of-view cutoff plate
- Compact, rugged, and low-cost
- Can be used as a surface reference for the QSP-Series underwater profiling instruments or PNF-300 Profiling Natural Fluorometer

Specifications

Optical Features Scalar Irradiance Collector:

1.9 cm (3/4") diameter solid Teflon® sphere optically connected to the main housing by a 2.5 cm aluminum-encased quartz light pipe.

Photodetector: Blue-enhanced high-stability silicon photovoltaic detector with dielectric and absorbing glass filter assembly.

PAR Spectral Response: Equal (better than ±10%) quantum response from 400 nm to 700 nm with response sharply attenuated above 700 nm and below 400 nm. Spectral response-induced errors will cause less than 5% errors in naturally occurring light fields.

Directional Response: Each instrument's directional response is optimized before final calibration. Front-to-side (approximately 85° from head on) response over all angles is equal (±6%) with response attenuated to 0 (at 95°) by the 20-cm diameter cutoff shield. Individual detector response plots are available as an option.

Sensitivity: When purchased alone, the sensor is calibrated in quanta/(cm²-sec))/volt. Nominal sensitivity is 1 volt = 1x10¹⁷ quanta/(cm²-sec) (slightly less than full sunlight). Noise level is typically less than 1 millivolt, temperature coefficient of the dark signal is less than 10 microvolts/°C, and response temperature coefficient is less than 0.15%/°C.

Electronic Features Measured Signals:

PAR Dynamic Range: 1.4x10⁻⁵ μE/(cm²·sec) to 0.5 μE/(cm²·sec)

Power Requirements:

Power: 12 VDC at 6 mA
Compatible Power Supplies:
DAS-186 automated Data
Acquisition System or QSP170BD Power Supply and Digital
Display (may be purchased separately).

Cable: QSC-240 shielded, weather-resistant cable with appropriate connectors (may be purchased separately).

Mechanical Features

Dimensions:

Diameter: 20.0 cm Height: 5.8 cm Weight: 0.9 kg

Housing Materials:

Collector: Hard anodized aluminum mounted with solid

PTFE sphere

Housing: Hard anodized alumi-

num

Environmental

Temperature Range: -2°C to 35°C

Calibration

The QSR-240 sensor is calibrated using a National Institute of Standards and Technology-(NIST) traceable 1000-watt type FEL Standard of Spectral Irradiance using procedures recommended by NIST. Annual recalibration is strongly recommended.

U.S. Patent No. 4,178,101



Biospherical Instruments Inc.

5340 Riley Street
San Diego, CA 92110-2621 USA
Phone: (619) 686-1888
Fax: (619) 686-1887
E-mail: sales@biospherical.com
URL: www.biospherical.com

Copyright © Biospherical Instruments Inc., 1999